

SUPPLEMENTARY REFERENCES

1. Group, C.C.P. The CCP4 Suite: Programs for Protein Crystallography. *Acta Crystallogr D* **50**, 760-763 (November 4, 1994).
2. Zamborlini, A. et al. Release of autoinhibition converts ESCRT-III components into potent inhibitors of HIV-1 budding. *Proc Natl Acad Sci U S A* **103**, 19140-5 (2006).
3. Shim, S., Kimpler, L.A. & Hanson, P.I. Structure/function analysis of four core ESCRT-III proteins reveals common regulatory role for extreme C-terminal domain. *Traffic* **8**, 1068-79 (2007).
4. Bajorek, M. et al. Biochemical Analyses of Human IST1 and Its Function in Cytokinesis. *Mol Biol Cell* (2009).
5. von Schwedler, U.K. et al. The protein network of HIV budding. *Cell* **114**, 701-13 (2003).
6. Babst, M., Wendland, B., Estepa, E.J. & Emr, S.D. The Vps4p AAA ATPase regulates membrane association of a Vps protein complex required for normal endosome function. *Embo J* **17**, 2982-93. (1998).
7. Yoshimori, T. et al. The mouse SKD1, a homologue of yeast Vps4p, is required for normal endosomal trafficking and morphology in mammalian cells. *Mol Biol Cell* **11**, 747-63. (2000).
8. Shim, S., Merrill, S.A. & Hanson, P.I. Novel Interactions of ESCRT-III with LIP5 and VPS4 and their Implications for ESCRT-III Disassembly. *Mol Biol Cell* **19**, 2661-72 (2008).
9. Studier, F.W. Protein production by auto-induction in high density shaking cultures. *Protein Expr Purif* **41**, 207-34 (2005).
10. Stuchell-Brereton, M.D. et al. ESCRT-III recognition by VPS4 ATPases. *Nature* **449**, 740-4 (2007).